

Diagnostic Slides (IFA)

Written by Administrator

Tuesday, 07 June 2011 22:14 - Last Updated Tuesday, 07 June 2011 23:52

Diagnostic Slides (IFA)

Teflon printed diagnostic slides for use in immunofluorescence test kits



Immuno-Cell is a leading manufacturer of Teflon printed diagnostic slides for use in IFA test kits. We have a long experience in this field and we manufacture these printed slides for the Diagnostic Industry all over the world.



Properties of the IFA slides

Diagnostic Slides (IFA)

Written by Administrator

Tuesday, 07 June 2011 22:14 - Last Updated Tuesday, 07 June 2011 23:52

For the production of our diagnostic slides for Immunofluorescence tests we use high quality microscopes

The hydrophobic Teflon paint is a in-house mix and is considered as a company secret.

The 2 main properties of these slides are :

a) the extreme hydrophobic Teflon printing which results in :

? no cross-contamination

type="margin: 0.0px 0.0px 0.0px 0.0px; font: 12.0px Arial;">? increased well capacity

b) the perfect wettable well surface which results in:

? a homogeneous spreading of the antigens within the well

? enhanced cell attachment on the well surface

Diagnostic Slides (IFA)

Written by Administrator

Tuesday, 07 June 2011 22:14 - Last Updated Tuesday, 07 June 2011 23:52

Our customers

a) Diagnostic Industry

For the manufacturers of test kits, most of the slides that we make are custom-made (OEM) following?

fications of our customers. These slides are printed with a design according to specifications, customized

a company logo, the name of the test, etc

b) Dealers of Diagnostic Products and Laboratory Supplies :

For our dealers we have an extended range of 80 standard executions. The most popular you find

and 5 of this catalogue.


These diagnostic slides are used in laboratories of Virology, Serology, Microbiology, Histology, etc

Diagnostic Slides (IFA)

Written by Administrator

Tuesday, 07 June 2011 22:14 - Last Updated Tuesday, 07 June 2011 23:52

The black printed standard executions are available in small quantities (1 box of 100 slides)



Teflon printed Slides

for Immunofluorescence tests

Our price policy for the IFA slides

Diagnostic Slides (IFA)

Written by Administrator

Tuesday, 07 June 2011 22:14 - Last Updated Tuesday, 07 June 2011 23:52

? Same prices for each design

? Prices are based on the quantity ordered per item number (not assorted).

? Same price for standard and custom-made executions, just the minimum order quantity changes

1) Standard executions (see page 4 and 5)

? Black Teflon printed minimum order : **100 slides**

? Other colors :. minimum order : **1.000 slides** (same price as standard executions)

available colors : blue - green - yeelow - grey - white - orange

? Orders without specification of the color will be printed automatically in black

2) Custom-made executions

Diagnostic Slides (IFA)

Written by Administrator

Tuesday, 07 June 2011 22:14 - Last Updated Tuesday, 07 June 2011 23:52

? Any design according to your specifications

? Color of your choice : black - blue - green - yellow - grey - white - orange - (or specific color on request)

minimum order : **2.000 slides** (same price as standard executions)

and a small contribution to the lay-out and screen costs (first order only)

The slides on which we print

? Our diagnostic slides are printed on the high quality "Pathology grade" slides with thickness 1
? 1.2mm (see

"MGF-Slides" page 16 in this catalogue)

Standard our diagnostic slides are printed on slides 76 x 26 mm (European size)

But on request we also print on 75 x 25mm slides (US size)

? The slides on which we print have 90° fine polished ground edges with "clipped corners"

Diagnostic Slides (IFA)

Written by Administrator


Tuesday, 07 June 2011 22:14 - Last Updated Tuesday, 07 June 2011 23:52

The "clipped corners" have the advantage that they do not damage the sealing foil when

they are vacuum packed after they have been coated by the manufacturer of the test kit.

? Depending of the printed design these slides are with or without a frosted end

Our available colors are:



Diagnostic Slides (IFA)

Written by Administrator

Tuesday, 07 June 2011 22:14 - Last Updated Tuesday, 07 June 2011 23:52

































Slide identification: the color will be printed automatically (e.g. black ? white ? red)

Diagnostic Slides (IFA)

Written by Administrator

Tuesday, 07 June 2011 22:14 - Last Updated Tuesday, 07 June 2011 23:52

 <p>61.100.12 36 wells 2 mm</p>	 <p>61.100.11 10 wells 5 mm</p>	 <p>61.100.17 10 wells</p>
 <p>61.100.33 12 wells 3 mm</p>	 <p>61.100.87 10 wells 5 mm</p>	 <p>61.100.29 12 wells</p>
 <p>61.100.50 30 wells 3 mm</p>	 <p>61.100.43 10 wells 5 mm HEp-2</p>	 <p>61.100.70 12 wells</p>
 <p>61.100.06 10 wells 4 mm</p>	 <p>61.100.77 12 wells 5 mm</p>	 <p>61.100.71 18 wells</p>
 <p>61.100.20 12 wells 4 mm</p>	 <p>61.100.26 14 wells 5 mm</p>	 <p>61.100.75 4 wells</p>
 <p>61.100.25 15 wells 4 mm</p>	 <p>61.100.36 18 wells 5 mm</p>	 <p>61.100.16 10 wells</p>
 <p>61.100.35 18 wells 4 mm</p>	 <p>61.100.60 4 wells 6 mm</p>	 <p>61.100.91 12 wells</p>
 <p>61.100.28 21 wells 4 mm</p>	 <p>61.100.05 8 wells 6 mm</p>	 <p>61.100.00 1 well</p>
 <p>61.100.32 24 wells 4 mm</p>	 <p>61.100.14 10 wells 6 mm</p>	 <p>61.100.47 1 well CHLAMYDIA</p>
 <p>61.100.01 6 wells</p>	 <p>61.100.02 10 wells</p>	 <p>61.100.03 2 wells</p>

Diagnostic Slides (IFA)

Written by Administrator

Tuesday, 07 June 2011 22:14 - Last Updated Tuesday, 07 June 2011 23:52

